# LuxaPrint Ortho

DMG Digital Enterprises SE  
Elbgaustraße 248 · 22547 Hamburg · Germany · www.dmg-dental.com  
093042/2020-09

## Instructions for use

English

### Product description

Light-curing resin based on (meth)acrylate for the production of dental devices. For DLP / SLA printers that work at wavelengths of 385 nm or 405 nm.

Devices with material parameters already defined in the system are recommended: e.g. DMG 3Demax/DMG 3Delite (DMG), D30/D20+/D20+ cartridge/D10+ (Rapid Shape), P20+/P10+ capsule (Straumann). Select material variants according to the wavelength used by the printer.

### Intended purpose

Manufacturing of drilling templates and orthodontic equipment using 3D printing.

### Limitations of use

Only insert LuxaPrint Ortho intra-orally in a fully polymerized state.

### Indications

- Surgical template: Tooth loss, in case of navigated implantology
- Orthodontic appliances: Stabilization of the tooth's position after the completion of active orthodontic treatment

### Contraindications

Do not use the material where there are known allergies to any of the components or contact allergies.

### Patient target group

Persons receiving treatment as part of a dental procedure.

### Intended users

Dentist, dental technician

### Notes for use

- During handling (e.g. transferring), the material should only be exposed to daylight / ambient light for as short a time as possible, since lengthy exposure to light can lead to unwanted curing of the material.
- Wear nitrile gloves, protective clothing and goggles when working with the material/printed object.

### For optimum properties of the end product:

- Ensure that the devices used during the production process are in proper condition and comply with machine and material parameters.
- Ensure clean working procedures. Residues on the devices can cause defects on the printed object.
- Use the material parameters already defined in the system for the cleaning units DMG 3Dewash (DMG), RS wash (Rapid Shape) and P wash (Straumann).
- Allow objects to cure completely.
- Use the material parameters already defined in the system for the light-curing units DMG 3Decure (DMG), RS cure (Rapid Shape) and P cure (Straumann).
- Minimum material thickness for structures: 1.5 mm.
- Maximum curing depth during post-curing: For solid objects and light-curing on both sides, the material can be up to 7 mm thick (at a curing depth of 3.5 mm).
- Thinner layers increase the accuracy of the printed object and increase printing time.
- Positioning the printed object as flat as possible on the printer's build platform. A flatter positioning increases accuracy and reduces printing time.

### Recommended use

**Note:** Observe the device and material parameters.

**Note:** Shake the bottle for at least 01:00 min before use.

**Note:** Maintain a processing temperature of 23 °C ± 2 °C (73.4 °F ± 3.6 °F).

**Note:** For optimal process reliability, scan the RFID code of the material (supported by DMG 3Demax/DMG 3Delite (DMG), D30/D20+/D20+ cartridge/D10+ (Rapid Shape); P20+/P10+ capsule (Straumann)).

### Finishing the printed object

Production steps	
Recommended drip-off time after raising the platform: 10:00 minutes	
<b>Note:</b> Finishing should take place immediately after the construction process where possible.	
1.	Carefully remove the printed object from the build platform. When using the printers DMG 3Delite (DMG), D10+ (Rapid Shape), P10+ capsule (Straumann), start directly from step 2.
2.	Cleaning: Clean the printed object with isopropanol (approx. 99 %) using an appropriate cleaning device and specified program (e.g. DMG 3Dewash, RS wash or P wash).
2. (Alternatively)	Pre-clean the printed object with ethanol (approx. 96 %) or isopropanol (approx. 99 %) in an ultrasonic cleaner for max. 03:00 minutes (if necessary with a brush). Optional: Clean openings, holes and gap areas with compressed air. Clean the printed object in a separate container with fresh ethanol (approx. 96 %) or isopropanol (approx. 99 %) for max. 02:00 minutes.
3.	Post-curing: Post-cure the printed object using the appropriate light-curing unit and the specified program (e.g. DMG 3Decure, RS cure or P cure).
4.	After removing the support, finish the surface of the printed object using the appropriate instruments and polishing paste.

### Preparation prior to use on the patient

#### Disinfection

According to the manufacturer's specifications, the following disinfectants may be used:

- PrintoSept-ID (on the basis of quart. ammonium salts)
- SprayActiv, alcoholic disinfectant spray (also contains didecyldimethyl-N-chloride)
- Dentavon (solution prepared from granulate; contains penta-potassium-bis (peroxymonosulfate)-bis(sulfate), anionic surfactants, non-ionic surfactants, soap, phosphonate)

#### Autoclavability

Before use on the patient, the printed object can be autoclaved once.

#### Autoclaving parameters for steam sterilization:

Temperature: 134 °C / 273 °F at a pressure of 2 bar and a duration of 05:00 min.

#### Interactions

Residues of methacrylate-based materials (uncured material residues) may influence the setting behavior of silicone impression materials.

#### Residual risks / side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

#### Warnings / precautions

- Keep out of the reach of children!
- Avoid skin contact with the liquid material and components prior to post-curing! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

## Performance characteristics of the device

Cured material	
Flexural strength	≥ 70 [MPa]
Transverse modulus	≥ 1 [GPa]
Shore-D hardness	≥ 60

Uncured material	
Viscosity	≤ 2.5 [Pa*s]
Light sensitivity to fluorescent lamps	≤ 05:00 [min]

## Storage / disposal

- Store in a dry place at room temperature (15°–25 °C/59°–77 °F) and protected from light!
- Even low exposure to light can trigger polymerization.
- Always keep the container tightly sealed; carefully seal immediately after every use.
- Do not use after the expiration date!
- The cured material can be discarded together with household waste following consultation with the disposal company and in compliance with regulations.
- Incinerate the resin mixture as hazardous waste in accordance with official regulations.

## Composition

EBPADMA, additives. Contains TPO.

**Note:** A glossary of all symbols displayed is at: <http://www.dmg-america.com/literature>

**Caution:** US Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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