DMG

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Instructions for use

English

Luxatemp-Automix Plus is a self-curing composite for fabricating temporary crowns and bridges, inlays, onlays and veneers. The two-component material based on multifunctional methacrylates is mixed automatically.

Indication

Temporary crowns and bridges, inlays, onlays, partial crowns, veneers and long-term temporaries

Timing

0:00 - 0:45 minutes

Insertion

2:00 - 3:00 minutes


Removal

6:00 - 7:00 minutes

End of setting

Automix system

1. Push up the small lever (1) at the back of the dispensing gun and pull back the slide (2) completely. Lift the plastic lever (3) on the top of the gun, insert cartridge and lock in by pushing down the lever (3). Ensure that the notch on the cartridge is aligned with the notch on the dispensing gun.
2. Remove the cartridge cap (4) by turning 90° anti-clockwise and discard it. Attach a new mixing cannula (5). Note the corresponding

notches. Turn the cannula 90° clockwise to lock it.

3. The material is dispensed through the mixing cannula by pulling the trigger (6) of the dispenser and can be applied directly.

When using a cartridge for the first time, express a pea-sized amount of material and discard.

Leave the used cannula on the cartridge to serve as a cap until the next use.

4. Push up the lever (1) at the back of the dispensing gun and pull the slide (2) back completely. Lift the plastic lever (3) and remove the cartridge.

Recommended use

1. Prior to preparation or extraction of a tooth, take a quadrant impression using alginate, silicone or a thermoplastic impression material.

2. Dry the prepared teeth and coat the stumps and surrounding tissue and any composite posts lightly with vaseline or a similar separating medium.
3. Select the desired Luxatemp shade and put the appropriate cartridge into the gun (see Automix system).
4. Dispense Luxatemp into your impression or plastic template. Apply the material to the occlusal surfaces first and then bring it gingivally, overbuilding it slightly. To prevent bubbles, it is important to keep the tip of the cannula immersed in the material.
5. Place the impression or plastic template filled with Luxatemp on the prepared teeth within 45 seconds with gentle pressure and hold firmly in place.
6. Remove the temporary from the mouth during the elastic phase, that is, 2:00 to 3:00 minutes after mixing.

The setting reaction has to be monitored intraorally (e.g. with a scaler) as

the temperature in the mouth has a significant influence on the setting time and the temporary can be removed only during the elastic phase.

7. The crown or bridge can be contoured and polished about 4:00 minutes after it is removed from the mouth. This is facilitated by removing the oxygen-inhibited layer on the surface of the temporary (e.g. with alcohol). The temporary can be shaped with slow-speed acrylic burs and discs etc. Use appropriate safety measures. Check occlusion and adjust if necessary with suitable instruments.
8. Varnish the surface with a light-curing high gloss varnish or polish manually with a rubber polisher or goat hair brush.
9. Cement the temporary using any temporary cement. Automatically mixing cements are recommended, e.g. DMG products such as TempoCem (zinc oxide/eugenol cement) or TempoCemNE (eugenol-free zinc oxide

cement), which can be applied directly into the temporary.

Cements containing eugenol may inhibit polymerization of resin-based luting cements.

10. Check occlusion.

Notes

- In molar areas with absent teeth, a ridge for pontic contact between the prepared teeth can be created by carving in the interproximal areas.
- When using a plastic template, the model is made from the impression. A plastic matrix is fabricated on this using a vacuum form machine.
- The precision of large temporaries can be improved further by placing them back in the impression after removing them from the mouth until setting is complete.

Repair of temporaries

When repairing a temporary, before applying Luxatemp, roughen the surface of the break point slightly with a bur or abrasion unit and apply a primer or e.g. Dry Coat according to the instructions for use. When repairing a recently fabricated temporary, it is not necessary to roughen it and use a primer.

1. Apply Luxatemp to all surfaces.
2. Join the parts and hold them together for 3:00 minutes.
3. After 6:00 minutes finish the repaired area as usual.

To speed up polymerization, the repaired temporary may be put in warm water at 50 °C for a few minutes. Alternatively, the Luxatemp temporary can also be repaired with the free-flow and light curing composite LuxaFlow. You will find information on handling this product in the LuxaFlow instructions for use.

Technical data

| | |
|---|------------|
| Compressive strength | 250 MPa |
| Transverse strength | 91.5 MPa* |
| Diametrical tensile strength | 36 MPa |
| Barcol hardness | 37 |
| Water absorption | 0.76 wt.%* |
| Maximum temperature during polymerization | ≈ 38 °C |

*Rzanny, A.; Welker D.; Göbel, R.: Phillip J., 13, 357-366 (1996)

Side effects

There have been no known systemic side effects to date. However, you should not use the material where allergies to one of the components are known to exist or in the rare cases where contact allergies exist.

Avoid contact with skin and eyes. In case of accidental contact, irrigate immediately with copious water and consult a doctor if necessary.

Composition

Glass filler in a matrix of multifunctional methacrylates; catalysts, stabilizers, additives. Free from methyl methacrylate and peroxides. Filler proportion:
44 wt.% = 24 vol.% (0.02-2.5 µm)

Storage

Store in a dry place below 25 °C. Use only at room temperature (refrigerated material is more viscous and takes longer to cure).

Do not use after the expiry date.

Presentations

Introductory pack

1 cartridge containing 76 g paste, colour A2

15 Automix-Tips

1 Automix-Dispenser

REF 110400

Refill pack

1 cartridge containing 76 g paste

| | |
|---------------------|------------|
| 15 Automix-Tips | |
| Colour A1 | REF 110401 |
| Colour A2 | REF 110402 |
| Colour A3,5 | REF 110403 |
| Colour B1 | REF 110404 |
| Colour Bleach Light | REF 110368 |

Economy pack

5 cartridges containing 76 g paste

| | |
|-----------------|------------|
| 75 Automix-Tips | |
| Colour A1 | REF 110405 |
| Colour A2 | REF 110406 |
| Colour A3,5 | REF 110407 |

Accessories

| | |
|----------------------------------|------------|
| 1 Automix-Dispenser Type 50 10:1 | REF 110411 |
| 1 Automix-Plunger Type 50 10:1 | REF 110412 |
| 45 Automix-Tips Blue 10:1 | REF 110409 |

Caution: Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in

which he or she practices to use or order the use of this device.

Keep away from children. For dental use only.

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