
# StatusBlue

---

MixStar

### PRODUKTBESCHREIBUNG

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Polysiloxane, Siliciumdioxid, Pigmente, Additive, Platinkatalysator

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### PRODUCT DESCRIPTION


StatusBlue-MixStar is an addition curing silicone impression material with neutral taste of medium consistence for automatic mixing and dispensing units (please refer to the respective instructions for use of the corresponding equipment). Its good recovery after deformation, good pourability and non slumping properties make StatusBlue-Mixstar an advantageous alternative to alginates.

### INDICATIONS

- Situation impressions of any kind
- Impressions for the preparation of provisionals
- Impressions for orthodontic models
- Impressions of the opposing jaw
- Impressions for the manufacture of model cast dentures


### PREPARATION

1


Slide locking mechanism into open position.  
**Do not remove the locking mechanism from the cartridge!**

2


Remove and discard the cap on the openings of base and catalyst outlets and attach a new mixing tip. Seat mixing tip **completely**.

3

Secure the mixing tip by sliding the locking mechanism into locked position.

4

Insert the cartridge into the mixing unit.

- Fresh cartridges are readied by discarding the first 3 cm of impression material that is extruded.
- Leave used tips on the cartridge after use. They serve as a cap.

### IMPRESSION TAKING

Due to its finely balanced consistency, StatusBlue-MixStar may be used in all trays. To get better adhesion to the tray an adhesive may be applied. Please refer to the respective instructions for use. The prepared tray is loaded with the impression material and should be introduced into the mouth no later than 1:15 minutes after start of mixing. The position must be kept in the mouth for at least 1:45 minutes. After removal, the impression is cleaned with water and disinfected.

### DISINFECTION AND CLEANING

Impressions can be disinfected with standard disinfection solutions.

### REMOVAL AND CASTING

After removal disinfect the impression. Impressions are then dried and stored at room temperature (max. 25 °C/77 °F). They can be poured 30 minutes after removal from the mouth. We recommend to use special purpose hard plaster of Paris or model cement.

### MANUFACTURE OF PROVISIONALS

Any provisional C&B can be used. The setting of these materials is somewhat faster than with alginates.

### PRECAUTIONS

- Do not mix with condensation or polyether type impression materials.
- Never allow the impression to come into contact with organic solvents or liquids, containing same, as this may result in swelling of the material and therefore may lead to inaccuracies.

### TECHNICAL DATA

Total working time (23 °C)	1:15 min
Recommended time in mouth]	1:45 min
Dimensional change after 24 h	≤ 0.4 %
Compression set	≤ 0.7 %
Strain in compression	≈ 7.0 %

### COMPOSITION

Addition curing polysiloxanes, silicon dioxide, food pigments, additives, platinum catalyst

### CLASSIFICATION

StatusBlue-MixStar fulfills the requirements of ISO 4823:2000, type 2: medium bodied.

### STORAGE

Do not store above 25 °C/ 77 °F.  
Do not use after expiration date.

### PACKAGING

1 Cartridge @ 390 ml paste, 10 MixStar-Tips	REF 909693
5 Cartridges @ 390 ml paste, 10 MixStar-Tips	REF 909694
50 MixStar-Tips	REF 903321

**Keep away from children!  
For dental use only!**

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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**DMG Chemisch-Pharmazeutische Fabrik GmbH**