



ClearShield 5% Sodium Fluoride Varnish

Instructions for Use

DESCRIPTION: ClearShield is a fluoride varnish which has a strong desensitizing action when applied to dental surfaces. ClearShield may be used on moist teeth and leaves a film of varnish which sets in the mouth and may allow visual control and verification.

INDICATIONS: ClearShield is a fluoride containing preparation for the treatment of dentin hypersensitivity, and for the reduction of post operative sensitivity. 1mL of ClearShield fluoride varnish contains: 50 mg of sodium fluoride, equivalent to 22.6 mg fluoride ion.

DOSAGE: For the treatment of dentin hypersensitivity, two applications at several days interval are sufficient in most cases. The fluoride content of ClearShield is dosed in such a way that neither acute nor chronic side effects are to be expected if applied according to instructions.

INSTRUCTIONS TO APPLY CLEARSHIELD

1. Peel open the foil lidding material
2. Remove brush and mix varnish until it reaches an even consistency
3. For best results, dry teeth prior to application by rubbing a 2x2 (or other disposable cotton sponge) along the tooth surfaces to remove any existing bio-film
4. Apply a very thin coat of varnish evenly to all of the desired tooth surfaces
5. Let the varnish dry for approximately 10 seconds and then instruct the patient to close their mouth
6. Varnish hardens on contact with saliva so the patient may leave immediately after application
7. When dry, any excess varnish may appear "spotty" and may slide off the tooth surface. This will not affect the efficacy of the product as there is still a full layer of varnish on the tooth
8. It is recommended that the patient be instructed to only eat soft foods and drink cold liquids for two hours after treatment

CAUTION: Federal law restricts the use of this device to sale by or on the order of a medical or dental professional. For use in the oral cavity.

CONTRAINDICATIONS: Ulcerative gingivitis and stomatitis.

INTERACTIONS: When ClearShield is applied, other fluoride preparations such as fluoride gels, should not be administered during the same day. The routine use of fluoride tablets should be interrupted for several days after treatment.

ADVERSE REACTIONS: Edematous swellings have been reported only in rare instances in some fluoride varnish products, especially after application to extensive surfaces. Dyspnea, although extremely rare, has occurred in asthmatic children. Nausea has been reported when extensive applications have been made. If required, varnish film can be removed with a thorough brushing.

SUPPLY: Individual packages containing 0.40 mL ClearShield 5% sodium fluoride varnish.

Reorder by calling 800-662-6383 or by fax 201-894-0213

ITEM#	For Reorder
799501	ClearShield 5% Sodium Fluoride Varnish Watermelon 35 pack
799502	ClearShield 5% Sodium Fluoride Varnish Watermelon 200 pack
799503	ClearShield 5% Sodium Fluoride Varnish Bubblegum 35 pack
799504	ClearShield 5% Sodium Fluoride Varnish Bubblegum 200 pack
799506	ClearShield 5% Sodium Fluoride Varnish Mint 35 pack
799507	ClearShield 5% Sodium Fluoride Varnish Mint 200 pack
799512	ClearShield 5% Sodium Fluoride Varnish Cookie Dough 35 pack
799513	ClearShield 5% Sodium Fluoride Varnish Cookie Dough 200 pack
799516	ClearShield 5% Sodium Fluoride Varnish Caramel 35 pack
799517	ClearShield 5% Sodium Fluoride Varnish Caramel 200 pack

Distributed by:
DMG America, LLC
242 South Dean Street
Englewood, NJ 07631
Phone #: 800-662-6383

MADE IN U.S.A.