# Ecosite

## Bulk Fill

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English

#### Product description

Ecosite Bulk Fill is a light-curing, radiopaque restorative material developed especially for use in the posterior region.

#### Indications

- Restorations in the posterior region (Classes I and II, including replacement of individual cusps)
- Restorations of class V (cervical caries, root erosions, wedge-shaped defects)
- Core build-up
- Restorations in deciduous teeth

#### Contraindications

- Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.
- Do not use directly on exposed pulp or in cases of allergies to any of the ingredients.

#### Interactions

- Avoid materials containing eugenol, moisture and oily air, as these may inhibit polymerization at the contact area.
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

#### Side effects

There are no known side effects to date. Hypersensitivity to components of the material cannot be ruled out in individual cases.

#### Basic safety information

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.

#### Notes for use

- ▶ For optimum polymerization, position the light-curing unit as close to the restorative material as possible and observe the recommendations relating to the light-curing time and the light intensity (see table). Observe the manufacturer's instructions for the light-curing unit used.
- ▶ Follow the manufacturer's instructions for other materials that are used with Ecosite Bulk Fill.
- ▶ If Ecosite Bulk Fill is applied into the patient's mouth straight from the Safetip, the Safetip must only be used on one patient for hygiene reasons.
- ▶ If a commercially available syringe gun is used to administer the contents of the Safetip composite tips, DMG advises following the manufacturer's instructions for the cleaning and disinfection of the syringe gun.

#### Recommended use

1. Clean the tooth. Remove all cleaning residue with water spray. It is advisable to use a rubber dam.
2. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Ensure that the areas close to the pulp are protected.
3. Condition the dental enamel and apply the adhesion primer according to the manufacturer's instructions for the material used. DMG advises using the universal adhesion primer Ecosite Bond.
4. If necessary, apply an underfilling of low-viscosity composite or compomer. We recommend using a matrix.

- Apply Ecosite Bulk Fill in the cavity with a maximum layer thickness of 5 mm for the contrast shade and 4 mm for the universal and light shades. Apply taller restorations one layer at a time. Light-cure each layer separately.

Light intensity	> 1,000 mW/cm <sup>2</sup>	> 600 mW/cm <sup>2</sup>
Curing time	20 s	40 s
Wavelength	420-480 nm	420-480 nm

**Note:** Additional buccal or lingual light-curing is advisable.

Remove metal matrices beforehand.

- Removal of excess material, contouring and fine processing are done using diamond finishing tools and flexible polishing discs of decreasing grain size, for example. Ecosite Bulk Fill can be given a high shine using polishing paste together with polishing cups or discs.

## Composition

Barium glass in a Bis-GMA-based matrix of dental resins.

Pigments, additives and catalyst. Filler content: 82% by weight = 65% by volume

The range of variation of the inorganic filler particles is between 0.02 and 0.7 µm

## Storage and shelf life

- Store at room temperature (2 to 25°C / 36 to 77°F).
- Close syringes / composite tip immediately after use so that the material does not cure in the ambient light.
- Do not use after the expiration date.

## Further information and specifications

The product specifications for Ecosite Bulk Fill material cured according to these instructions for use include the following values: a compressive strength of ≥ 300 MPa, a flexural strength of ≥ 120 MPa and a depth of cure of ≥ 4 mm.

**Caution:** U.S Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

**Note:** A glossary of all symbols displayed is at: <http://www.dmg-america.com/literature>

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