LuxaPrint Tray

DMG Digital Enterprises SE
Elbgaustraße 248 · 22547 Hamburg · Germany · www.dmg-dental.com
093041/2020-09

Instructions for use

English

Product description

Light-curing resin based on (meth)acrylate for the production of dental devices. For DLP / SLA printers that work at wavelengths of 385 nm or 405 nm. Select material variants according to the wavelength used by the printer.

Devices with material parameters already defined in the system are recommended: e.g. DMG 3Demax/DMG 3Delite (DMG), D30/D20+/D20+ cartridge/D10+ capsule (Rapid Shape), P20+/P10+ capsule (Straumann).

Intended purpose

Manufacturing of custom impression trays using 3D printing.

Limitations of use

Only insert LuxaPrint Tray intraorally in a fully polymerized state.

Indications

All dental situations that require final impressions, for instance, for the fabrication of fixed or removable prosthetics.

Contraindications

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended users

Dentist, dental technician, dental assistant staff

Notes for use

- During handling (e.g. transferring), the material should only be exposed to daylight / ambient light for as short a time as possible, since lengthy exposure to light can lead to unwanted curing of the material.
- Wear nitrile gloves, protective clothing and goggles when working with the material/printed object.

For optimum properties of the end product:

- Ensure that the devices used during the production process are in proper condition and comply with machine and material parameters.
- Ensure clean working procedures. Residues on the devices can cause defects on the printed object.
- Use the material parameters already defined in the system for the cleaning units DMG 3Dewash (DMG), RS wash (Rapid Shape) and P wash (Straumann).
- Allow objects to cure completely.
- Use the material parameters already defined in the system for the light-curing units DMG 3Decure (DMG), RS cure (Rapid Shape) and P cure (Straumann).
- Maximum curing depth during post-curing: For solid objects and light-curing on both sides, the material can be up to 7 mm thick (at a curing depth of 3.5 mm).
- Thinner layers increase the accuracy of the printed object and increase printing time.
- Position the printed object as flat as possible on the printer's build platform. A flatter positioning increases accuracy and reduces printing time.

Recommended use

Note: Observe the device and material parameters.

Note: Shake the bottle for at least 01:00 min before use.

Note: Maintain a processing temperature of 23 °C ± 2 °C (73.4 °F ± 3.6 °F).

Note: For optimal process reliability, scan the RFID code of the material (supported by DMG 3Demax/DMG 3Delite (DMG), D30/D20+/D20+ cartridge/D10+ capsule (Rapid Shape); P20+/P10+ capsule (Straumann)).

Finishing the printed object

Production step	
Recommended drip-off time after raising the platform: 10:00 minutes	
Note: Finishing should take place immediately after the construction process where possible.	
1.	Carefully remove the printed object from the build platform. When using the printers DMG 3Delite (DMG), D10+ (Rapid Shape), P10+ capsule (Straumann), start directly from step 2.
2.	Cleaning: Clean the printed object with isopropanol (approx. 99 %) using an appropriate cleaning device and specified program (e.g. DMG 3Dewash, RS wash or P wash).
2. (Alternatively)	Pre-clean the printed object with ethanol (approx. 96 %) or isopropanol (approx. 99 %) in an ultrasonic cleaner for max. 03:00 minutes (if necessary with a brush). Optional: Clean openings, holes and gap areas with compressed air. Clean the printed object in a separate container with fresh ethanol (approx. 96 %) or isopropanol (approx. 99 %) for max. 02:00 minutes.
3.	Post-curing: Post-cure the printed object using the appropriate light-curing unit and the specified program (e.g. DMG 3Decure, RS cure or P cure).
4.	After removing the support, finish the surface.

Disinfection / Sterilization

Note: Do not use heat-based disinfection or sterilization methods. Otherwise, the workpiece may be deformed.

Interactions

Residues of methacrylate-based materials (uncured residues) may influence the setting behavior of silicone impression materials.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings / precautions

- Keep out of the reach of children! Even low exposure to light can trigger polymerization.
- Avoid skin contact with the liquid material and components prior to post-curing! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Avoid breathing in the fumes.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Performance characteristics of the device

Cured material	
Flexural strength	≥ 70 [MPa]
Transverse modulus	≥ 1 [GPa]
Shore-D hardness	≥ 50

Uncured material	
Viscosity	≤ 2.5 [Pa*s]
Light sensitivity to fluorescent lamps	≤ 05:00 [min]

Storage / disposal

- Store in a dry place at room temperature (15°–25 °C/59°–77 °F) and protected from light!
- Even low exposure to light can trigger polymerization.
- Always keep the container tightly sealed; carefully seal immediately after every use.
- Do not use after the expiration date!
- The cured material can be discarded together with household waste following consultation with the disposal company and in compliance with regulations.
- Incinerate the resin mixture as hazardous waste in accordance with official regulations.

Composition

UDMA, TEDMA, additives.

Note: A glossary of all symbols displayed is at: <http://www.dmg-america.com/literature>

Caution: US Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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