LuxaFlow

DMG Chemisch-Pharmazeutische Fabrik GmbH
Elbgaustraße 248 · 22547 Hamburg · Germany
www.dmg-dental.com
091754/#2_2015-02

CE 0482

EN Instructions for use

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Instructions for use

English

Product description

LuxaFlow is a light-curing, free-flowing composite that is especially suited for the repair of Luxatemp or other bis-acrylate materials, as well as for fillings for cavity classes III, IV and V.

Indications

- Individualization and repair of temporary bis-acrylates
- Small fillings for cavity classes III, IV and V
- Minimally invasive fillings
- Underfillings

Contraindications

- Do not use the material in the event of allergies against one of the ingredients or in the event of contact allergies.
- Do not use the material if using a drain or the recommended application technique is not possible.
- Do not apply the material to the exposed pulp.

Basic safety instructions

- Only for use by dental professionals!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of unintended contact with the skin, wash the affected area thoroughly and immediately with soap and water.
- Avoid contact with the eyes! In the event of unintended contact with the eyes, flush thoroughly and immediately with plenty of water and consult a doctor if necessary.
- When working on temporary restorations, wear a mask and protective goggles; otherwise, grinding dust could enter the eyes and respiratory tract!

Indications for use

- Light-curing units should emit at 450 nm and be inspected regularly.
- The light intensity should total at least 400 mW/cm².
- Place the light-curing unit as close as possible to the material.
- Do not remove the oxygen-inhibited layer since it will be needed for bonding with the next layer.

- When using LuxaFlow in combination with other materials, adhere to the respective manufacturer's instructions as closely as possible.

Recommended use

Use as individualization and repair material

Note: During individualization or repair of a recently-produced temporary restoration, it is not necessary to abrade or use a bonding agent (proceed with step 2).

Caution! When working on temporary restorations, grinding dust could enter the eyes and respiratory tract!

▶ Wear a mask and protective goggles.

1. Repair of a temporary restoration that is already in place: Lightly abrade the pertinent spot on the restoration with a cutter or sandblaster. Apply a bonding agent in accordance with the directions for use.
2. Select the desired LuxaFlow shade. Remove the syringe cap and screw on the Luer-Lock cannula.
3. Press the syringe and apply the material in 2 mm layers (maximum) with the assistance of the Luer-Lock tip.
4. Light-cure each layer for 30 to 40 sec.
5. Shape and finish each area as customary.

Use as filling material

1. Clean and dry the preparation.
2. Prepare cavity in accordance with the general rules of the adhesive technique.

Caution! If the pulp is not adequately protected, there is a possible risk of damage to the pulp and pain for the patient.

▶ Protect the pulp in deep cavities.

3. Conduct etching in accordance with the preferred technique (e.g. Total Etch). Adhere to directions for use of the applied material!
4. Apply bonding agent (e.g. Contax, TECO). Observe corresponding directions for use!
5. Select the desired LuxaFlow shade. Remove the syringe cap and screw on the Luer-Lock cannula.
6. Press the syringe and apply the material in 2 mm layers (maximum) with the assistance of the Luer-Lock tip.
7. Light-cure each layer for 30 to 40 sec.
8. Remove excess material, contour and refine the filling surfaces, e.g. with a diamond bur and flexible discs of decreasing grain sizes.
9. Polish the filling surface and adjoining surface with polishing paste. Use polishing cups or polishing discs to polish.

Side effects

No systemic side effects are known to date. In isolated cases, hypersensitivity or contact allergy to components of the material cannot be excluded.

Interactions

Humidity, oil-contaminated air, and materials containing eugenol may prevent polymerization at the contact point.

▶ Avoid humidity, oil-contaminated air, and materials containing eugenol!

Discoloration may occur with the use of cationic mouthwashes, plaque revealers or chlorhexidine.

Composition

Dental glass in an optimized matrix made of BIS-GMA. Filler content 63 wt.% = 41 vol.%. The range of the inorganic filler particles is between 0.02 and 3 µm.

Classification

ISO 4049:1988; depth of cure (40 seconds): 3.5 mm.

Storage and shelf life

- ▶ Store in a dry place at room temperature (15 to 25 °C/59 to 77 °F)!
- ▶ After use, reseal syringes!
- ▶ Do not use past the expiry date!

Exclusively distributed in the USA and Canada by: DMG America LLC Englewood, New Jersey 07631, U.S., Phone: (800) 662-6383.

Emergency: Infotrac (800) 535-5053

Caution: US Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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