PermaCem 2.0

Self-Adhesive

DMG Chemisch-Pharmazeutische Fabrik GmbH
Elbgaustraße 248 · 22547 Hamburg · Germany
www.dmg-dental.com
092757/#1_2015-01

EN Instructions for use

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Instructions for Use

English

Product description

PermaCem 2.0 is a universal, self-adhesive and dual-cure composite luting cement for the permanent cementation of crowns, bridges, inlays, onlays and root posts. The use of an additional enamel-dentine bonding agent or of a metal-zirconia primer is not required.

Indications

- Permanent cementation of inlays, onlays, crowns, bridges and root posts composed of metal and metal alloys (precious and non-precious metals), metal ceramics and composites, all-ceramics such as zirconium oxide, aluminum oxide and all forms of silicate ceramics (e.g. lithium disilicate, feldspar), as well as fiber-reinforced composites (root posts)
- Permanent cementation of crowns and bridges on implant abutments comprised of titanium or zirconia

Contraindications

Do not use the material:

- If the preparation does not offer sufficient retention (e.g. as in the case of veneers)
- In cases of allergies to any of the ingredients or in cases of contact allergies
- For direct capping of the pulp tissue

Basic safety instructions

- ▶ For dental use only!
- ▶ Keep out of the reach of children!
- ▶ Avoid contact with the skin! In the event of inadvertent skin contact, immediately wash the affected area thoroughly with soap and water.
- ▶ Avoid contact with the eyes! In the event of inadvertent contact with the eyes, immediately rinse the eyes thoroughly with plenty of water and consult a physician if necessary.
- ▶ In the event of application to dentine in close proximity to pulp, use pulp protection.

Side effects

There are no known side effects to date. Hypersensitivity to components of the material cannot be excluded in individual cases.

Interactions

Products containing eugenol, for example temporary cements or temporary root canal filling materials, may impair curing of the composite luting cement.

- ▶ Do not use products containing eugenol and other substances that inhibit polymerization together with the material. If possible, use eugenol-free cements (e.g. TempoCem NE).
- ▶ Make sure that any residual cement is completely removed from the preparation or core following removal of the temporary restoration.
- ▶ From the beginning of pretreatment until permanent cementation, ensure that neither the restoration nor the tooth surfaces become contaminated.

Timing

≈1:00 minute	Working time (including mixing time)*
0:30–1:00 minutes after insertion of the restoration/root post	Removal of excess (without light curing)
1–2 seconds	Pre-light curing time to make removal of excess more feasible (optional)
7:00 minutes	Maximum intraoral curing time (including mixing time)
At least 20 seconds	Light-curing time (optional)

*The time stated applies at a room temperature of 23°C. Higher temperatures may shorten these times, while lower temperatures may increase them.

Instructions for use

- ▶ Ambient light may lead to premature curing of the composite luting cement. You must therefore protect PermaCem 2.0 from ambient light after expressing it from the Smartmix syringe.
- ▶ Ensure that excess cement in hard-to-reach areas (proximal, marginal gingival areas) is removed promptly.
- ▶ The PermaCem 2.0 curing process can be accelerated. In order to do so, cure PermaCem 2.0 using a suitable curing light.
- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm². Place the light as close as possible to the material.

Recommended use

For luting indirect restorations such as crowns, bridges, inlays and onlays

1. Prepare the clean and dry restoration as follows:
 - ▶ Etchable silicate ceramic: etch the restoration with 5 % hydrofluoric acid in accordance with the instructions from the manufacturer of the etching gel and rinse thoroughly using a water spray. Subsequently dry with air that is free of oil and water, and prepare the surface using a silanizing agent (e.g. Vitique Silane) in accordance with the manufacturer's instructions.
 - ▶ Oxide ceramics, metals and metal alloys: pretreat in accordance with the manufacturer's instructions. Unless otherwise stated, sandblast the inside of the restoration with aluminum oxide, clean with alcohol, and dry with oil- and water-free air.
 - ▶ Other materials: pretreat in accordance with the manufacturer's instructions.
2. If required, rinse the cavity or preparation thoroughly using a water spray to remove any residual temporary cement, and then dry. In doing so, ensure that the dentine retains a little residual moisture and does not become overly dry.
3. Apply PermaCem 2.0 inside the pretreated restoration.
4. Insert the restoration within 1:00 minute after mixing has begun and attach it to the prepared tooth using slight pressure.
5. Remove excess cement within 0:30 to 1:00 minute of inserting the restoration with the aid of a microbrush, brush, foam pellet or scaler. Remove interdental excess carefully using dental floss.
It is alternatively possible to remove excess cement by briefly curing it using light (1–2 seconds) and then removing it with a scaler.
6. Instruct the patient to adopt habitual occlusion.
7. Allow the material to chemically cure fully for 7:00 min from the start of mixing.
8. In the case of translucent restorations, additionally cure the material using a suitable curing light for at least 20 seconds.

Cementation of root posts

1. Endodontically prepare the root canal, clean it, and dry it using a paper point. Ensure that the dentine retains a little residual moisture and does not become overly dry.
2. Prepare the clean and dry root post as follows:


- ▶ Fiber-reinforced root post (e.g. LuxaPost): Clean the root post with ethanol, dry with air that is free of oil and water, and apply silane (e.g. Vitique Silane) where appropriate in accordance with the manufacturer's instructions.
 - ▶ Root posts from other materials: pretreat in accordance with the manufacturer's instructions.
3. Apply PermaCem 2.0 directly into the prepared canal with the aid of the appropriate tip. In doing so, insert the tip as deeply as possible into the canal.

Note: During the application of the composite luting cement, ensure that the tip remains submerged in the material and that it is only moved upwards with the rising material.


4. Insert the root post within 1:00 minute after mixing has begun and cement using slight pressure.
5. Remove the excess cement within 0:30 seconds to 1:00 minute of inserting the root post with the aid of a microbrush, brush, foam pellet or scaler.
6. Allow the material to chemically cure fully for 7:00 min from the start of mixing.

Note: The material can be additionally cured using a suitable curing light for at least 20 seconds.

Using the Smartmix syringe


1. Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90° counterclockwise, and discard it.

2. Attach a new mixing tip [B].

Note: In order to achieve an optimum mixing result, DMG recommends using the mixing tips available from DMG. All mixing tips are also available in refill packs.

Note: Make sure that the notches on the mixing tip line up with those on the Smartmix syringe [C].

3. Turn the mixing tip clockwise 90° to lock it.
4. The material is mixed in the mixing tip when it is expressed and can then be applied directly.

Note: When using a syringe for the first time, squeeze out a pea-sized amount of material and discard it!

After application, leave the used mixing tip on the Smartmix syringe. It serves as a cap!

Composition

Barium glass in a Bis-GMA-based matrix of dental resins, pigments, additives and catalysts. Filler content: 69 % by weight = 51 by vol. % (0.02 – 3.0 μm).

Classification

ISO 4049:2009 (including color stability).

Storage and shelf life

- ▶ Store in a dry place at temperatures of 2 to 8°C (36 to 46°F)!

- ▶ Protect against exposure to direct sunlight!
- ▶ Do not use after the expiration date!

Caution: U.S. Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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