# LuxaPrint Ortho Plus

**DMG Digital Enterprises SE**

Elbgaustraße 248 · 22547 Hamburg · Germany · www.dmg-dental.com  
093046/2021-09

## Instructions for use

English

### Product description

Light-curing resin based on (meth)acrylate for the production of dental devices. For DLP / SLA printers that work at wavelengths of 385 nm or 405 nm.

Devices with material parameters already defined in the system are recommended: e.g. DMG 3Demax/DMG 3Delite (DMG), D30/D20+/D20+ cartridge/D10+ (Rapid Shape), P20+/P10+ capsule (Straumann). Select material variants according to the wavelength used by the printer.

### Indications for use

Fabrication of orthodontic and dental appliances such as mouthguards, nightguards, splints and positioners.

### Limitations of use

- Only insert LuxaPrint Ortho Plus intraorally in a fully polymerized state.
- LuxaPrint Ortho Plus is not suitable for fabrication of sport mouthguards.

### Contraindications

Do not use the material where there are known allergies to any of the components or contact allergies.

### Patient target group

Persons receiving treatment as part of a dental procedure. Orthodontic and dental appliances such as mouthguards, nightguards, splints and positioners are used after active orthodontic treatment and/or to protect the teeth and reduce damage caused by bruxing or nighttime grinding and prevention of the noise associated with bruxing, grinding and snoring.

### Intended users

Dentist, dental technician

### Notes for use

- During handling (e.g. transferring), the material should only be exposed to daylight / ambient light for as short a time as possible, since lengthy exposure to light can lead to unwanted curing of the material.
- Polymerized LuxaPrint Ortho Plus material cannot be reused. Unpolymerized material can be stored in the printer tank or filtered (recommended mesh size 190 µm) and returned to the bottle.
- Wear nitrile gloves, protective clothing and goggles when working with the material/printed object.

### For optimum properties of the end product:

- Ensure that the devices used during the production process are in proper condition and comply with machine and material parameters.
- Ensure clean working procedures. Residues on the devices can cause defects on the printed object.
- Use the material parameters already defined in the system for the cleaning units DMG 3Dewash (DMG), RS wash (Rapid Shape) and P wash (Straumann).
- Allow objects to cure completely.
- Use the material parameters already defined in the system for the light-curing units DMG 3Decure (DMG), RS cure (Rapid Shape) and P cure (Straumann).
- Minimum material thickness for structures: 1.5 mm.
- Maximum curing depth during post-curing: For solid objects and light-curing on both sides, the material can be up to 7 mm thick (at a curing depth of 3.5 mm).
- Thinner layers increase the accuracy of the printed object and increase printing time.
- Position the printed object as flat as possible on the printer's build platform. A flatter positioning increases accuracy and reduces printing time.

### Recommended use

**Note:** Observe the machine and material parameters. Do not use LuxaPrint Ortho Plus with any device, resin materials, scanners, printers, post-curing units, CAD/CAM software, templates, and tools, other than those specifically identified in the LuxaPrint Ortho Plus labeling.

**Note:** Shake the bottle for at least 01:00 min before use.

**Note:** Maintain a processing temperature of 23 °C ± 2 °C (73.4 °F ± 3.6 °F).

**Note:** For optimal process reliability, scan the RFID code of the material (supported by DMG 3Demax/DMG 3Delite (DMG), D30/D20+/D20+ cartridge/D10+ (Rapid Shape); P20+/P10+ capsule (Straumann)).

### Finishing the printed object

Production steps	
Recommended drip-off time after raising the platform: 10:00 minutes	
<b>Note:</b> Finishing should take place immediately after the construction process where possible.	
1.	Carefully remove the printed object from the build platform. When using the printers DMG 3Delite (DMG), D10+ (Rapid Shape), P10+ capsule (Straumann), start directly from step 2.
2.	Cleaning: Clean the printed object with isopropanol (approx. 99 %) using an appropriate cleaning device and specified program (e.g. DMG 3Dewash, RS wash or P wash).
2. (Alternatively)	Pre-clean the printed object with ethanol (approx. 96 %) or isopropanol (approx. 99 %) in an ultrasonic cleaner for max. 03:00 minutes (if necessary with a brush). Optional: Clean openings, holes and gap areas with compressed air. Clean the printed object in a separate container with fresh ethanol (approx. 96 %) or isopropanol (approx. 99 %) for max. 02:00 minutes.
3.	Post-curing: Post-cure the printed object using the appropriate light-curing unit and the specified program (e.g. DMG 3Decure, RS cure or P cure).
4.	After removing the support, finish the surface of the printed object using the appropriate instruments and polishing paste.

### Interactions

Residues of methacrylate-based materials (uncured material residues) may influence the setting behavior of silicone impression materials.

### Residual risks / side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out

### Warnings / precautions

- Keep out of the reach of children!
- Avoid skin contact with the liquid material and components prior to post-curing! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Avoid breathing in the fumes.
- Follow the manufacturer's instructions for other products that are used with the material/materials.

- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.
- Do not use heat-based disinfection or sterilization methods. Otherwise, the workpiece may be deformed.

## Performance characteristics of the device

Cured material	
Flexural strength	≥ 70 [MPa]
Transverse modulus	≥ 1 [GPa]
Shore-D hardness	≥ 60

Uncured material	
Viscosity	≤ 2.5 [Pa*s]
Light sensitivity to fluorescent lamps	≤ 05:00 [min]

## Storage / disposal

- Store in a dry place at room temperature (15°–25 °C/59°–77 °F) and protected from light!
- Even low exposure to light can trigger polymerization.
- Always keep the container tightly sealed; carefully seal immediately after every use.
- Do not use after the expiration date!
- The cured material can be discarded together with household waste following consultation with the disposal company and in compliance with regulations.
- Incinerate the resin mixture as hazardous waste in accordance with official regulations.

## Composition

EBPADMA, additives. Contains TPO.

**Note:** A glossary of all symbols displayed is at: <http://www.dmg-america.com/literature>

**Caution:** US Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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