


Icon Smooth Surface

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English

Product description

Icon Smooth Surface is an innovative product for the micro-invasive treatment of lesions. It was especially developed for smooth surface infiltration and is particularly suitable for the treatment and masking of white spots.

Indications for use

- Micro-invasive treatment of non-cavitated cariogenic lesions.
- Esthetic improvement of non-cavitated lesions caused by demineralization and/or mild to moderate fluorosis (on smooth surfaces)

Contraindications

Do not use the material where there are known allergies to any of the components or contact allergies.

Limitations of use

Sufficient drying and protection of gingiva is not possible

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended users

Dentist

Notes for use

- Light-curing units should emit at 450 nm and should be checked regularly.
- The light intensity should be a minimum of 800 mW/cm².
- Place the light unit as close to the material as possible.
- The etched surface must not be touched or contaminated with saliva or moisture before the treatment resumes. If contamination occurs after drying, re-etch for approx. 10 seconds and repeat drying with Icon Dry.
- For effective treatment results a sufficiently dry working area is critical. Therefore, appropriate measures for complete isolation must be taken. Do not use rubber dams made from thermoplastic elastomers such as Flexi-Dam (Coltène Whaledent/Hygenic).
- Discard syringes after use.
- The syringes included in the treatment unit contain sufficient material for treating two to three smooth surface lesions.
- Multiple smooth surface lesions can be treated at the same time.
- For topical white spots an area of 2 mm beyond the lesion site is etched. At the discretion of the dentist, the entire smooth surface area can be etched and infiltrated in the case of extensive white spots occurring after bracket removal.
- The overall esthetic result may be improved by means of a whitening treatment before infiltration, especially in the case of fluorosis treatments. The manufacturer's instructions must be taken into account.

Recommended use (for image, see »Handling«)

1. Preparation

- Before the start of treatment, clean the affected tooth and adjacent teeth. Remove all cleaning residue with water spray.
- Ensure that the working area is sufficiently dry. For this purpose, appropriate measures such as a rubber dam or liquid dam should be taken. Make sure that the lesion to be treated remains fully accessible.

2. Exposing the lesion body

- Pretreatment measures: To ensure successful infiltration, the complete hypermineralized surface layer must be removed prior to applying Icon Infiltrant in order to gain access. Similarly, access to deeper white spots should also be gained first.
- The use of suitable assistive devices (e.g., medium grit abrasive disk, microabrasion or sandblasting) is optional prior to etching.

Application of Icon Smooth Surface Etch:

- Screw the Vestibular Tip onto the Icon Smooth Surface Etch syringe.
- Apply an ample amount of Icon Smooth Surface Etch onto the lesion site by turning the syringe shaft carefully, let it take effect for 02:00 minutes and activate the effect by moving the Tip occasionally. Remove excess material with cotton roll (Fig. 1).
- Aspirate off Icon Smooth Surface Etch and rinse with water for at least 30 seconds. Then dry using oil-free and water-free air.

3. Result preview and drying

For effective treatment results it is necessary to dry the lesion to its full depth.

In this step a preview of the final result is shown.

Application of Icon Dry:

- Screw the Application Tip onto the Icon Dry syringe, apply an ample amount of material onto the lesion and allow it to take effect for 30 seconds. Dry carefully with oil-free and water-free air (Fig. 2).
- When wetted with Icon Dry, the whitish-opaque coloration on the etched enamel should fade. If this is not the case, repeat step 2.
- After a visual check, dry the lesion with oil-free and water-free air to remove the Icon Dry.

4. Infiltration

Do not apply Icon Infiltrant under direct operating light as this may cause the material to set prematurely.

Application of Icon Infiltrant

- Screw a new Vestibular Tip onto the Icon Infiltrant syringe.
- Apply an ample amount of Icon Infiltrant by turning the shaft.
- Allow Icon Infiltrant to infiltrate for 03:00 minutes, occasionally activate the infiltration by moving the syringe and top up if necessary (Fig. 3).
 In case of deeper and larger defects the esthetic result can be improved by extending the exposure time.

- Remove excess material with a cotton roll and, if necessary, dental floss.
- Light-cure Icon Infiltrant for a minimum of 40 seconds from all sides.
- Repeat the application with Icon Infiltrant as previously described but with a contact time of 01:00 minute
- Remove the rubber dam.
- Use polishing cups, or similar, for the surface finish.

Interactions

Do not use together with substances which inhibit polymerization. This is especially valid for materials which contain eugenol, such as temporary cements.

Residual risks / side effects

- In some cases, contact allergies with similar composite products have been reported.
- Contact of Icon Smooth Surface Etch with the oral mucosa causes a white surface coloration. This will subside after a few days.

Warnings / precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Caution! Icon Smooth Surface Etch can cause chemical burns if not used properly

- Avoid contact with soft oral tissue, eyes, and skin. In case of inadvertent contact rinse with copious amounts of water, and, if necessary, consult a physician.
- Limit the application to the surface to be treated. Avoid contact of the etching gel with gingiva, exposed dentin, root element and adjacent teeth.
- Use only the enclosed syringe and application tip to apply the etching gel in the patient's mouth.
- Discard the applicator tip after use. Replace syringe cap.

Performance characteristics of the device

Penetration coefficient: approx. 147 [cm/s]

Storage / disposal

- Store at 2-25 °C (36–77 °F) in a dry place.
- Protect against exposure to direct sunlight!
- Do not use after the expiration date!

Composition

- Icon Smooth Surface Etch: Hydrochloric acid (approx. 47 %), water (approx. 33 %), glycerol (approx. 15 %), pyrogenic silicic acid, surface-active substances
- Icon Dry: Ethanol (approx. 99 %)
- Icon Infiltrant: Methacrylate-based resin matrix (approx. 99%), initiators, additives

Note: A glossary of all symbols is displayed at: <http://www.dmg-america.com/literature>

Caution: US Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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