

00:01:00

ONE MINUTE PRODUCT MEETING

Vitique

Vitique is a premium resin cement developed especially to enhance the production of perfect veneers. It offers several innovations unmatched by other veneer cements, such as superior color stability for long lasting beautiful veneers, consistent ability to reflect the shade of the veneer and facilitate the natural transparency of the tooth, and optimum flow characteristics allowing for uncomplicated placement of the thinnest veneer.

Vitique's lack of sensitivity to ambient light allows for relaxed, stress-free placement of multiple veneers.

FEATURES & BENEFITS

- › Perfect color coordination
- › Try-In + cement ensures a precise preview of final case
- › Specially designed ribbon Veneer tip for convenience
- › Remove excess material with ease

KIT TYPES

- › The Cementation Kit: cements and Try-In pastes
- › The Accessory Kit: etching gel, silane and DMG's unique bonding agent, TECO Total-Etch



PROCEDURES

- › Inlays › Onlays › Veneers

COMPANION PRODUCTS

- › Luxatemp

ITEM # & DESCRIPTION

213080		Veneer System
213081		Cementation System
213082		Accessory Kit
213038		A1 Cement
213039		A2.5 Cement
213040		A4 Cement
213041		B1 Cement
213042		BL Cement
213043		Transparent Cement
213044		White Cement
213045		Pink Cement
213048		A1 Try-in Paste
213049		A2.5 Try-in Paste
213050		A4 Try-in Paste
213051		B1 Try-in Paste
213052		BL Try-in Paste
213053		Transparent Try-in Paste
213054		White Try-in Paste
213055		Pink Try-in Paste

LOOKING FOR MORE INFO?

Visit www.dmg-america.com/products.