

00:01:00

ONE MINUTE PRODUCT MEETING



DentaMile Desk MC system

A fusion of simplicity, practicality, and intuitive design, the new system, DentaMile Desk MC, delivers highly-precise and impressive performance, in a compact size.



DentaMile Desk MC-5

The DentaMile Desk MC-5 DLP printer fits anywhere, yet offers a spacious build platform for producing multiple devices. Print accurate mouthguards, splints, nightguards, models and bleaching trays effortlessly thanks to the high precision 2K light source.

DentaMile Wash MC

The DentaMile Wash MC offers simple, preset programs that achieve fast, automatic cleaning of 3D-printed devices. The generous container, magnetic mixer, and removable basket ensure isopropyl alcohol effectively cleans prints and the unit.

DentaMile Cure MC

Equipped with dual-frequency Omni 360 degree technology, the DentaMile Desk Cure MC ensures uniform curing through two wavelengths for a variety of printing resins.

Features & Benefits

- › User-friendly system
- › DentaMile connect ready
- › Highly precise printed devices
- › Indicated for a wide range of 3D dental applications
- › Validated workflow
- › Compatible with all standard dental CAD programs (STL files)
- › Compact design
- › Unprecedented quality
- › Increased production
- › Reproducible precision
- › Coordinated processes
- › Intuitive operability
- › Predictable outcomes
- › Guaranteed regulatory compliance

Companion products

- › DentaMile connect
- › LuxaPrint resins

Ordering Information

- 70500 | DentaMile Desk MC-5
- 170501 | DentaMile Wash MC
- 170502 | DentaMile Cure MC
- 170510 | DentaMile Desk MC-5 Resin Reservoir
- 170511 | DentaMile Desk MC-5 Magnetic Platform
- 170512 | DentaMile Desk MC-5 Flex Plate

For more product information visit:

dmg-america.com/3dprinting

1-800-662-6383 • info@dmg-america.com

A smile ahead together

