

# 00:01:00

## ONE MINUTE PRODUCT MEETING

### LuxaPrint

LuxaPrint products are high-quality, light-cured resins for use with digital light processing printers that work at a wavelength of 385nm. Each material combines fast printing with outstanding accuracy for all your 3D printing needs.

#### LuxaPrint Tray

Indicated for the fabrication of custom impression trays.

#### LuxaPrint Model

Indicated for the fabrication of dental models that are suitable for the production of thermoformed splints.

#### LuxaPrint Ortho

Indicated for the fabrication of custom surgical guides and splints with the highest precision.

#### LuxaPrint Ortho Plus

Indicated for the manufacturing of orthodontic and dental appliances such as mouthguards, nightguards, splints and tooth positioners.

#### LuxaPrint Ortho Flex

A biocompatible material for the additive production of flexible splints like bleaching trays or indirect bonding trays.

#### LuxaPrint Cast

Indicated for the fabrication of residue-free burn-out molds (e.g. model and denture framework castings, crowns and bridges).

#### LuxaPrint Gingiva

Ideal for creating gingival masks to help simulate the periodontal status for dental laboratory models or master models.

### FEATURES & BENEFITS

- › High quality
- › Light cured
- › Easy finishing and processing
- › Compatible with DMG 3Dmax printer
- › Validated for use with Straumann printers
- › Wide range of applications

### ITEM # & DESCRIPTION

171264	LuxaPrint Model	1-200gm Bottle, Beige (385/405nm)
171294	LuxaPrint Model	1-200gm Bottle, Grey (385/405nm)
171200	LuxaPrint Model	1-500gm Bottle, Beige (385/405nm)
171292	LuxaPrint Model	1-500gm Bottle, Grey (385/405nm)
171201	LuxaPrint Model	1-1,000gm Bottle, Beige (385/405nm)
171293	LuxaPrint Model	1-1,000gm Bottle, Grey (385/405nm)
171266	LuxaPrint Tray	1-200gm Bottle, Turquoise (385/405nm)
171204	LuxaPrint Tray	1-500gm Bottle, Turquoise (385/405nm)
171205	LuxaPrint Tray	1-1,000gm Bottle, Turquoise (385/405nm)
171262	LuxaPrint Cast	1-200gm Bottle, Dark Red (385/405nm)
171202	LuxaPrint Cast	1-500gm Bottle, Dark Red (385/405nm)
171203	LuxaPrint Cast	1-1,000gm Bottle, Dark Red (385/405nm)
171256	LuxaPrint Ortho	1-200gm Bottle, Transparent (385nm)
171206	LuxaPrint Ortho	1-500gm Bottle, Transparent (385nm)
171207	LuxaPrint Ortho	1-1,000gm Bottle, Transparent (385nm)
171259	LuxaPrint Ortho Plus	1-200gm Bottle, Transparent (385nm)
171210	LuxaPrint Ortho Plus	1-500gm Bottle, Transparent (385nm)
171211	LuxaPrint Ortho Plus	1-1,000gm Bottle, Transparent (385nm)
171273	LuxaPrint Ortho Flex	1-200gm Bottle, Transparent (385nm)
171250	LuxaPrint Ortho Flex	1-500gm Bottle, Transparent (385nm)
171252	LuxaPrint Ortho Flex	1-1,000gm Bottle, Transparent (385nm)
171278	LuxaPrint Gingiva	1-200gm Bottle, Pink (385nm)
171240	LuxaPrint Gingiva	1-500gm Bottle, Pink (385nm)
171241	LuxaPrint Gingiva	1-1,000gm Bottle, Pink (385nm)



### LOOKING FOR MORE INFO?

Visit [www.dmg-america.com/3dprinting](http://www.dmg-america.com/3dprinting)