

# LuxaBite®

## Innovative rigid bisacryl bite registration material.

- Combines superior handling properties and a firmer set to produce the most accurate and stable bite registration available.

LuxaBite's rigid bisacryl chemistry withstands temperature variations and handling better than any other material, so the bite registration you take is the same one the lab receives! Unlike traditional plaster and polyvinyl siloxane bite registration materials, LuxaBite exhibits no compression or flexing when mounted, ensuring total predictability. In fact, LuxaBite is a valuable tool for checking the accuracy of plaster models and you can even use it to mount several sets of casts to the same bite.

### Ordering Information

Item #	Description
110560	LuxaBite Refill Kit (1-50ml Cartridge, 15 Automix Tips)
110609	LuxaBite Mixing Tips (Bag of 25)
110411	Type 50 4:1/10:1 Applicator Gun (1 Gun)

LuxaBite Refill Automix Kit 110560



# O-Bite®

## Rigid-setting bite registration material with proprietary vinyl polysiloxane formulation.

- Shore A Hardness of 94 provides unrivaled stability and resistance to breakage.

O-Bite is a superior new bite registration material that provides exceptional accuracy, predictability and speed. Specifically developed to yield more stable and precise bite registrations, O-Bite's proprietary vinyl polysiloxane formulation offers a Shore A Hardness of 94 providing unrivaled stability and resistance to breakage.

O-Bite makes bite registrations easier for you, your patients and your lab. O-Bite is a soft, non-gritty, thixotropic, patient-friendly orange material with a neutral taste and pleasant orange scent. O-Bite's Automix Type 50 dispenser system allows for easy mixing and dispensing, precise placement and maximum efficiency, with Automix tips that creates 32% less material waste. O-Bite requires no special technique and will not stick to undercuts or existing restorations. Once you've removed the registration from the mouth, simply trim it with a standard scalpel, disinfect it and send it off to the lab.

### Ordering Information

Item #	Description
999765	O-Bite Refill Kit (4-50ml Cartridge, 15 Automix Tips)
999202	O-Bite Mixing Tips (Bag of 50)
999507	Type 50 1:1 Applicator Gun (1 Gun)

O-Bite Automix kit 999765

