


TECO™ Total Etch Control

Total etch. Totally easy.

■ The one push-button bonding system for all light-cured materials.



Advantages

- Outstanding bonding properties.
- Unique fast and easy application system.
- Proven Total Etch technology.

 **DMG**
AMERICA



Hold TECO between thumb and index finger.



Push the silver button.



TECO after activation.



Remove the Microbrush with the ready-to-use bonding material...



... and apply easily to enamel and dentin.

TECO™ Total Etch Control.

Reliable bonding with push-button convenience.

Get the reliability you need with the dependability of a total etch bonding system in an innovative single-dose dispenser with TECO.

Outstanding bonding properties.

TECO offers a microtensile bond strength of 38 Mpa* – far superior to the accepted industry standard for optimal adhesion. Its non-acetone formulation works with all curing lights.

Unique fast and easy application system.

The bonding material is activated by a simple push on the patented “SilvR” single-dose package and may be used immediately. The built-in micro brush applicator easily slides out of the SilvR dose capsule. Application is a one-step, one-hand, fast, clean and safe process.

Proven Total Etch technology.

TECO can be used with a wet or moist bonding technique. The Total Etch method ensures exceptional bond strength and durability.

*University of Alabama study, 2007

Ordering Information: TECO™ Total Etch Control.

213320 60 TECO unit-dose (0.12 mL) packages of Bond and 1 syringe (2 mL) Etching Gel – Medium Viscosity



For more information, call us at **1-800-662-6383**. To find out more about all of our innovative dental products, visit our website at www.dmg-america.com

DMG America

242 South Dean Street
Englewood, New Jersey 07631
P: 800-662-6383 • 201-894-5505

F: 201-894-0213 • www.dmg-america.com

