



## 4 Ways LuxaBite's Unsurpassed Rigidity Delivers Unsurpassed Bite Registrations

### 1. Exceptional stability and handling

LuxaBite is automatically mixed to a homogenous, bubble-free, non-slumping consistency and can be directly applied, which prevents the material from flowing into approximal areas.

### 2. Resistance to compression and flexing

LuxaBite's high hardness (Barcol 25 after 1 hour) prevents compression or flexing when mounted.

### 3. Precise reproductions for correction-free restorations

LuxaBite provides the prerequisite for a precise reproduction of the oral situation, enabling the dentist to integrate perfectly fitting restorations without the need for corrections.

### 4. Impossible to break but easy to trim

The contrasting blue color record is break-proof and easy to trim.



### The triumph of rigidity

**LuxaBite** combines superior handling properties and a firmer setting to yield the most accurate, hard and stable bite registrations on the market. Due to its bisacryl chemistry, it exhibits better mechanical properties than VPS bite registration materials.