



5 Ways StatusBlue Outperforms Alginates

1. Dimensional stability

Impressions taken with StatusBlue maintain their dimensions; they neither shrink nor swell.

2. Pourability at any time

StatusBlue impressions don't have to be kept moist, so there's no "Pour it or lose it!" pressure. You can pour the model today, tomorrow or next week! You can even cold sterilize StatusBlue impressions and re-pour as needed.

3. No time-consuming hand-mixing

StatusBlue puts an end to the mess, trouble and lost chair time caused by having to hand-mix alginates.

4. Easy removal

StatusBlue allows quick, distortion-free removal from dental devices.

5. Slump-free tray dispensing

StatusBlue has an innovative A-silicone formula for a firm consistency and better patient acceptance than alginates.



The state-of-the-art alginate replacement material

StatusBlue's innovative A-silicone formula makes it quick and easy to take situation impressions for temporaries, model-cast dentures, orthodontic models, and opposing jaw models. **Its time-saving prep and setting times will save nearly 60 hours of valuable chair time for the typical dentist every year compared to alginates.**