

# 00:01:00

## ONE MINUTE PRODUCT MEETING

### Constic

#### Etch, bond and fill, in one step!

Constic is a light-cured, self-etching, self-adhesive, radiopaque, flowable restorative composite indicated for a diverse set of clinical indications.

This 3-in-1 flowable provides the dentist with a clinically flexible, multi-purpose tool for various applications; Class I and Class II restorations, an adhesive liner under composite restorations, as a pit and fissure sealant, or for repairs and modifications to existing restorations or provisional restorations and many more.



### FEATURES & BENEFITS

- > Self-etching and self-adhesive
- > Faster, easier, more efficient treatment process
- > Radiopacity and tooth-like fluorescence
- > Prevents postoperative sensitivities
- > High bond strength
- > Low elastic modulus (high flowability and less stress on the tooth)

### PROCEDURES


- > Small restorations of class I and small occlusal primary tooth cavities
- > Pit and fissure sealing
- > Base lining of class I and class II restorations
- > Blocking out and filling of undercuts

### BOND STRENGTH

- > Enamel: 20.23 Mpa
- > Dentin: 24.23 Mpa

### ITEM # & DESCRIPTION

220700   A1	220704   B1
220701   A2	220705   Opaque White
220702   A3	220712   Accessories Pack
220703   A3.5	



### LOOKING FOR MORE INFO?

Visit [www.dmg-america.com/products](http://www.dmg-america.com/products).