

00:01:00

ONE MINUTE PRODUCT MEETING

Constic

Etch, bond and fill, in one step!

Constic is a light-cured, self-etching, self-adhesive, radiopaque, flowable restorative composite indicated for a diverse set of clinical indications.

This 3-in-1 flowable provides the dentist with a clinically flexible, multi-purpose tool for various applications; small Class I and Class III restorations, an adhesive liner under composite restorations, as a pit and fissure sealant, or for repairs and modifications to existing restorations or provisional restorations and many more.



FEATURES & BENEFITS

- > Self-etching and self-adhesive
- > Faster, easier, more efficient treatment process
- > Radiopacity and tooth-like fluorescence
- > Prevents postoperative sensitivities
- > High bond strength
- > Low elastic modulus (high flowability and less stress on the tooth)

PROCEDURES


- > Small restorations of class I and class III, small occlusal primary tooth cavities
- > Pit and fissure sealing, extended fissure sealing
- > Cavity lining of class I and class II restorations
- > Blocking out and filling of undercuts

BOND STRENGTH

- > Enamel: 20.23 Mpa
- > Dentin: 24.23 Mpa

ITEM # & DESCRIPTION

220700 A1	220704 B1
220701 A2	220705 Opaque White
220702 A3	220712 Accessories Pack
220703 A3.5	



LOOKING FOR MORE INFO?

Visit www.dmg-america.com/products.