

# 00:01:00

## ONE MINUTE PRODUCT MEETING

### MiniDam

DMG's MiniDam is a latex-free, silicone-based gingival protection device that protects the proximal area during treatment – quickly, easily and most importantly, comfortably for the patient.

Its innovative design ensures that the site being treated is kept relatively dry, while protecting the gingiva from the materials that could cause harm, such as etching gel. The MiniDam can be placed in just a few seconds and stabilizes itself without the use of clamps, while still allowing open access to the treatment area.



#### FEATURES & BENEFITS

- › Reliable local protection of the proximal region
- › Quick, easy and comfortable for patients
- › Secure hold without clamps
- › Latex-free
- › Ideal for bleaching, acid etching or Icon resin infiltration

#### PROCEDURES

- › Composite restorations
- › Resin infiltration

#### COMPANION PRODUCTS

- › Icon Proximal

#### ITEM # & DESCRIPTION

220381 | MiniDam (20 pieces)



#### LOOKING FOR MORE INFO?

Visit [www.dmg-america.com/products](http://www.dmg-america.com/products).