

00:01:00

ONE MINUTE PRODUCT MEETING

LuxaPrint

LuxaPrint products are high-quality, light-cured resins for use with digital light processing printers that work at a wavelength of 385nm. Each material combines fast printing with outstanding accuracy for all your 3D printing needs.

LuxaPrint Tray

Indicated for the fabrication of custom impression trays.

LuxaPrint Model

Indicated for the fabrication of dental models that are suitable for the production of thermoformed splints.

LuxaPrint Ortho

Indicated for the fabrication of custom surgical guides and splints with the highest precision.

LuxaPrint Cast

Indicated for the fabrication of residue-free burn-out molds (e.g. model and denture framework castings, crowns and bridges).

LuxaPrint Gingiva

Ideal for creating gingival masks to help simulate the periodontal status for dental laboratory models or master models.

LOOKING FOR MORE INFO?

Visit www.dmg-america.com/3dprinting

FEATURES & BENEFITS

- › High quality
- › Light cured
- › Easy finishing and processing
- › Compatible with DMG 3Demax printer
- › Validated for use with Straumann printers
- › Wide range of applications

ITEM # & DESCRIPTION

171264	LuxaPrint Model	1-200gm Bottle, Beige (385/405nm)
171294	LuxaPrint Model	1-200gm Bottle, Grey (385/405nm)
171200	LuxaPrint Model	1-500gm Bottle, Beige (385/405nm)
171292	LuxaPrint Model	1-500gm Bottle, Grey (385/405nm)
171201	LuxaPrint Model	1-1,000gm Bottle, Beige (385/405nm)
171293	LuxaPrint Model	1-1,000gm Bottle, Grey (385/405nm)
171266	LuxaPrint Tray	1-200gm Bottle, Turquoise (385/405nm)
171204	LuxaPrint Tray	1-500gm Bottle, Turquoise (385/405nm)
171205	LuxaPrint Tray	1-1,000gm Bottle, Turquoise (385/405nm)
171262	LuxaPrint Cast	1-200gm Bottle, Dark Red (385/405nm)
171202	LuxaPrint Cast	1-500gm Bottle, Dark Red (385/405nm)
171203	LuxaPrint Cast	1-1,000gm Bottle, Dark Red (385/405nm)
171256	LuxaPrint Ortho	1-200gm Bottle, Transparent (385nm)
171206	LuxaPrint Ortho	1-500gm Bottle, Transparent (385nm)
171207	LuxaPrint Ortho	1-1,000gm Bottle, Transparent (385nm)
171278	LuxaPrint Gingiva	1-200gm Bottle, Pink (385nm)
171240	LuxaPrint Gingiva	1-500gm Bottle, Pink (385nm)
171241	LuxaPrint Gingiva	1-1,000gm Bottle, Pink (385nm)

