



Constic

Do more with less.



Constic: Do More With Less

- Self-etching, self-adhesive flowable composite
- Combines an etching gel, bonding agent and flowable composite in one for a faster, easier procedure







Features & Benefits

- ↪ Self-etching and self-adhesive
- ↪ Faster, easier and more efficient treatment process
- ↪ Radiopacity and tooth-like fluorescence
- ↪ Prevents postoperative sensitivities
- ↪ Minimizes potential source of mistakes
- ↪ High bond strength
- ↪ Low elastic modulus





Self Adhesion, Secure Bonding, Simple Handling



*Source: Muñoz-Viveros CA (2012); Data on file, DMG, Hamburg, Germany.



Functional Shade System





Indications

- ↗ Small restorations of class I and small occlusal primary tooth cavities
- ↗ Pit and fissure sealing
- ↗ Base lining of class I and class II restorations
- ↗ Blocking out and filling of undercuts



Constic is applied to the cavities/fissures using the Luer-Lock-Tip



Constic is massaged into the entire tooth surface for 25 s using the included brush



**Simple application,
great impact. Constic.**

Apply directly, quickly
and precisely.



Product Information

DMG America Item #	Product & Contents
220700	Shade A1 (2 Syringes @ 2g Composite, 1 Brush Holder, 20 Clip-on Brushes, 20 Luer-lock Silver Tips)
220701	Shade A2 (2 Syringes @ 2g Composite, 1 Brush Holder, 20 Clip-on Brushes, 20 Luer-lock Silver Tips)
220702	Shade A3 (2 Syringes @ 2g Composite, 1 Brush Holder, 20 Clip-on Brushes, 20 Luer-lock Silver Tips)
220703	Shade A3.5 (2 Syringes @ 2g Composite, 1 Brush Holder, 20 Clip-on Brushes, 20 Luer-lock Silver Tips)
220704	Shade B1(2 Syringes @ 2g Composite, 1 Brush Holder, 20 Clip-on Brushes, 20 Luer-lock Silver Tips)
220705	Shade Opaque White (2 Syringes @ 2g Composite, 1 Brush Holder, 20 Clip-on Brushes, 20 Luer-lock Silver Tips)

Product Images

