



Honigum Pro

Producing True Impressionist Masterpieces



Impressions with superior accuracy



From the makers of Luxatemp and LuxaCore, Honigum Pro is the only VPS that goes where you want it, stays where you place it and has the unique ability to thin under load. It makes its way into tiny crevices and hard-to-reach areas that other impression materials can miss.

Honigum Pro's true thixotropic formula does not slump or run, but remains stable where you place it. It also exhibits remarkable strength to resist tearing upon removal, preventing the need for retakes. Because, the better the impression, the better the restoration.

Plus, DMG's "Snap-Set" technology provides an unbeatable combination of generous working time, even while still in the mouth, while offering a conveniently short intraoral time.



FEATURES & BENEFITS:


- Exhibits dynamic viscosity when placed under pressure
- Flexible working time
- Rapid snap-set technology
- Maximum tear strength
- Extremely hydrophilic
- Subtle honey aroma

APPLICATIONS:

- Crowns
- Bridges
- Inlay/Only
- Full-mouth restorations
- Implants

Stays Where You Place it.


Viscosity at rest*
(under no pressure) Low shear - 10 Pa



At low shear stress Honigum Pro Light remains stable, which means no running from the preparation.

Goes Where You Need it.

Viscosity under load*
(pressure upon insertion) High shear - 1000 Pa



At high shear stress Honigum Pro Light exhibits superior flow properties to record the finest details.

"I love Honigum Pro for the accuracy and reliability it gives me. I rarely have to take a second impression with it. The margins are super clear. For me, this is everything I want in an impression material."
– Sarah Winter, DMD
Private Practice, La Jolla, CA

*Source : Internal DMG measurements, data on file.



Stability and flowability have met their match

With an ideal combination of flowability and stability, flowing well around implant impression posts and transfer caps, Honigum Pro remains at the application site. It is the perfect material for soft-tissue, crown and bridge as well as inlay and onlay impression taking.

Perfect Precision For All Types of Impressions

Honigum Pro Light is a light consistency wash material with extraordinary stability, improved flowability, high toughness and precise detail.

Honigum Pro Heavy is a heavy consistency tray material that exhibits extreme accuracy, stability, easy removal and rigidity for locked-in detail.

Save Time with Faster-Setting Quad Fast

Honigum Pro Light Quad Fast has perfect stability, improved flowability and flawless detail in a light consistency, **faster-setting material**.

Honigum Pro Heavy Quad Fast has extreme accuracy, stability, easy removal and rigidity for locked-in detail. This heavy consistency, **faster-setting material** is ideal for quadrant and single unit impressions.



Gingival formers after the implants have healed into position.



Syringing around the impression posts.



Complete syringing around impression posts.



Fixation of impression posts in an individual, open impression tray, provides positionally stable fixation.

Excellent Shaping & Recovery After Deformation

Honigum Pro-Putty offers optimal shaping by hand with its putty consistency, has **outstanding recovery** after deformation, excellent dimensional stability and tear strength.

Precise Impressions in a Single Step

Honigum Pro Mono is a medium consistency monophasic impression material that is precise, easy to remove, has optimal flowability and high tear strength, all in a **single step**.



“Honigum Pro’s handling characteristics and results meet all of my requirements for an exceptional impression material.”

– Gary DeWood, DDS
SPEAR Education, Scottsdale, AZ



The MixStar eMotion: Clean, easy, and hands-free

Efficient and dynamic for mixing and dispensing pre-filled 380ml cartridges, the MixStar eMotion includes two built-in timers that can indicate working and setting times with up to 99 different combinations.

The Advanced Blue MixStar Mixing Tip has a redesigned cannula for a better mix with 11% less material waste.

Ordering Information

Honigum Pro

Item#	Product	Contents
989831	Honigum Pro Light	2-50ml Cartridges, 10 Automix and Intra-Oral Tips
989539	Honigum Pro Light	4-25ml Cartridges, 10 Automix & Intra-Oral Tips
989832	Honigum Pro Light Quad Fast	2-50ml Cartridges, 10 Automix and Intra-Oral Tips
989639	Honigum Pro Light Quad Fast	4-25ml Cartridges, 10 Automix & Intra-Oral Tips
989766	Honigum Pro Heavy	4-50ml Cartridges, 8 Automix Tips
989537	Honigum Pro Heavy	1-380ml Cartridge, 10 MixStar Mixing Tips
989838	Honigum Pro Heavy Quad Fast	4-50ml Cartridges, 8 Automix Tips
989638	Honigum Pro Heavy Quad Fast	1-380ml Cartridge, 10 MixStar Tips
989830	Honigum Pro Mono	4-50ml Cartridges, 8 Automix Tips & Intra-Oral Tips
989568	Honigum Pro Mono	1-380ml Cartridge, 10 MixStar Tips
989374	Honigum Pro Putty	1-380ml Cartridge, 10 MixStar Tips

Tips and Accessories

Item#	Product	Contents
999482	MixStar eMotion Machine	1-Mixing Machine Unit
993334	MixStar Blue Mixing Tips	Box of 50 Tips
999400	Honigum Pro Automix Tips	Bag of 50 Tips (Heavy, Mono)
110424	Pro Automix/Intra-Oral Mixing Tips	Bag of 35 Each (Light)
999490	Intra-Oral Tips	Bag of 50 Tips (Mono)
999394	Tray Adhesive	1-10ml Bottle
110253	Type 25 1:1 Applicator Gun	1-Gun
999507	Type 50 1:1 Applicator Gun	1-Gun



DMG America, LLC
65 Challenger Road, Suite 340, Ridgefield Park 07660
tel: 800.662.6383 email: info@dmg-america.com
web: www.dmg-america.com

