



# LuxaBite

Innovative rigid bisacryl bite registration



# The triumph of rigidity

LuxaBite combines superior handling properties and a firmer setting to yield the most accurate, hard and stable bite registrations on the market. Due to its bisacryl chemistry, it exhibits better mechanical properties than VPS bite registration materials.

## Stability

LuxaBite is automatically mixed to a homogenous, bubble-free, non-slumping consistency and can be directly applied. This prevents the material from flowing into approximal areas. The high hardness (Barcol 25 after 1h) of LuxaBite prevents compression or flexing when mounted.

## Precise Reproduction

With its outstanding rigidity, LuxaBite provides the prerequisite for a precise reproduction of the oral situation enabling the dentist to integrate perfectly fitting restorations without the need for corrections. The contrasting blue color record is break-proof and easy to trim.

## Ordering Information

Item #	Product	Contents
110560	LuxaBite Refill Kit	1-50ml Cartridge, 15 Automix Tips
110609	LuxaBite Mixing Tips	Bag of 25
110411	Type 50 4:1/10:1 Applicator Gun	1-Gun



LuxaBite applied to a row of teeth (max. 3 to 4 teeth).



Patient in occlusion.



Open and close bite situation after 45 to 60 seconds.



Easy excess removal with crossed-toothed burr.

DMG America, LLC  
65 Challenger Road, Suite 340, Ridgefield Park, NJ 07660  
tel: 800.662.6383 email: info@dmg-america.com  
web: www.dmg-america.com

