



# O-Bite

Exceptional accuracy, predictability  
and speed



# Superior bite registration

O-Bite is a superior bite registration material that provides you with exceptional accuracy, predictability and speed. Specifically developed to yield more stable and precise bite registrations, O-Bite's proprietary vinyl polysiloxane formulation offers a Shore A Hardness of 94 providing unrivaled stability and resistance to breakage.

## High final hardness

O-Bite has been designed to be the hardest VPS-based bite registration material. This property eliminates errors caused by compression of softer bite materials during articulation on the model. The superior hardness and perfect flow of O-Bite insures exact replication of intercuspal position for perfectly fitting restorations.

## Excellent break resistance

O-Bite exhibits outstanding break resistance. The structure of the VPS material is achieved through silicone chains of equal length which create a stable linking formation when cured. This feature guarantees high break resistance during removal, cutting and repositioning on the model.

## Ordering Information

Item #	Product	Contents
999765	O-Bite	4-50ml Cartridges, 15 Automix Tips
999202	O-Bite Automix Tips	Bag of 50
999507	Type 50 1:1 Applicator Gun	1 Gun



Direct application of O-Bite.



Patient in desired occlusion.



Bite registration before removal.



Easy to cut with a scalpel or a silicone cutter.

DMG America, LLC  
65 Challenger Road, Suite 340, Ridgefield Park, NJ 07660  
tel: 800.662.6383 email: info@dmg-america.com  
web: www.dmg-america.com

