

00:01:00

ONE MINUTE PRODUCT MEETING

LuxaBite

LuxaBite is an innovative rigid bite registration material that combines superior handling properties and a firmer setting, to yield the most accurate, hard and stable bite registration on the market. Its bisacryl chemistry yields better mechanical properties than other VPS bite registration materials. It is automatically mixed to a homogeneous, bubble-free, non slumping consistency that can be directly applied.

The stability of LuxaBite prevents the material from flowing into approximal areas, and its superior final hardness prevents compression or flexing when mounted. Plus, the contrasting blue color record is break-proof and easy to trim.



LOOKING FOR MORE INFO?
Visit www.dmg-america.com/products.

tel: 800.662.6383

email: info@dmg-america.com

web: www.dmg-america.com

FEATURES & BENEFITS

- › The hardest bite registration material on the market (Barcol hardness 25 after 1 hour)
- › No compression or flexing when mounted
- › Firm set for superior accuracy
- › Excellent stability
- › Easy to trim with a high-speed bur
- › Break-proof
- › Minimal resistance to closure of teeth
- › Taste and odor free
- › Contrasting color: Blue

PROCEDURES

- › Bite registration

COMPANION PRODUCTS

- › Luxatemp

ITEM # & DESCRIPTION

110560 | LuxaBite

110609 | Automix Tips

