



YOUR ALLY IN AMALGAM ALLOY



The alloy type **you want.**
The consistent mix and
optimum handling **you need.**

YOUR ALLY IN AMALGAM ALLOY

The DMG alloy product line has the alloy type you want that meets all your amalgam alloy needs.



Premium Dispersed Phase Amalgam Alloy (Regular & Fast Set)

The Premium line ensures good marginal integrity and eliminates open contacts while providing longer lasting, corrosion resistant restorations.

Item numbers

	Single Spill (400mg)		Double Spill (600 mg)		Triple Spill (800mg)	
	Regular	Fast Set	Regular	Fast Set	Regular	Fast Set
50 Capsule Jar	#61008	#61009	#61012	#61013	#61016	#61017
500 Capsule Jug	#61020	#61021	#61024	#61025	#61028	#61029
	Green/Grey	Green/White	Dark Blue/Grey	Dark Blue/White	Yellow/Grey	Yellow/White



Composition Contains: 50% Mercury by Weight

Approximate Alloy Content: Silver 53%, Tin 23%, Copper 24%, Zinc<.5%

FEATURES AND BENEFITS:

- Ensures good marginal integrity
- Eliminates open contacts
- Provides longer lasting and corrosion resistant restorations
- Requires no pinching or twisting of the capsule to open
- Self-activating capsule



Type-T Spherical Amalgam Alloy

The Type-T Spherical Amalgam Alloy is the strongest type of alloy made that resists fracture and ditching. Its fast setting allows for immediate packing and carving.

Item numbers

	Single Spill (400mg)	Double Spill (600 mg)	Triple Spill (800mg)
	50 Capsule Jar	#61039	#61040
500 Capsule Jug	#61035	#61036	#61037
	Aqua/Royal Blue	Aqua/Gray	Aqua/White



Composition Contains: 42.5% Mercury by Weight

Approximate Alloy Content: Silver 56%, Tin 29%, Copper 15%

FEATURES AND BENEFITS:

- Strongest type of alloy made
- Fast setting
- Fracture resistant
- Immediate silky carve
- Self-activating capsule
- Color-coded capsules



Manufactured by: DMG America LLC • 242 South Dean Street • Englewood, NJ 07631
Phone: (800) 662-6383 • www.dmg-america.com