



DentaMile Desk MC system

Embrace digital dentistry effortlessly

A 3D printing fusion of **simplicity, practicality, and intuitive** design

Jump into digital dentistry with the validated DentaMile Desk MC 3D printing system. The printer, wash, and cure units deliver highly-precise dental appliances, in a perfectly compact size.

Start printing today with the DentaMile Desk MC-5 printer, DentaMile Wash MC, and DentaMile Cure MC post-processing units.

FEATURES AND BENEFITS

- ✔ User-friendly printer, wash and cure units
- ✔ Accurately printed devices
- ✔ Perfect for nightguards, models, bleaching trays and more
- ✔ DentaMile connect ready
- ✔ Validated workflow
- ✔ Compatible with all standard dental CAD programs (STL files)
- ✔ Compact design that fits anywhere



DENTAMILE DESK MC-5

With its compact dimensions, the DentaMile Desk MC-5 DLP printer fits anywhere, yet offers a spacious build platform for producing multiple devices. Print accurate mouthguards, splints, nightguards, models and bleaching trays effortlessly thanks to the high precision 2K light source. Plus, it's DentaMile connect ready.

TECHNICAL DATA

Size (width/depth/height)	11.81 x 15.74 x 16.63 in. / 35.27 lbs
Work platform	5.11 x 2.87 x 3.54 in.
Technology	DLP @ 385 nm
Native pixel resolution	68 µm
Connections	2 USB ports, 1 LAN port 1 main connection outside
Additional information	<ul style="list-style-type: none"> • 2K light engine • 5" color touch display • WLAN and LAN • High-precision image alignment and orientation correction process



DENTAMILE CONNECT

Networking people and devices is the goal of DentaMile connect. Experience the convenience of cloud-based design software where step-by-step, guided workflows make it possible to efficiently integrate digital dentistry into your daily routine.



LUXAPRINT

A perfect companion to the DentaMile Desk MC system, the LuxaPrint resin portfolio includes a suite of light-cured materials for high-quality dental device fabrication with the utmost precision.



DENTAMILE WASH MC

The DentaMile Wash MC offers a simple wash program to achieve fast, automatic cleaning of your 3D-printed devices. The generous container, magnetic mixer, and removable basket ensure isopropyl alcohol effectively cleans your devices.

TECHNICAL DATA

Size (width/depth/height)	12.59 x 12.99 x 9.68 in.
Capacity	4.67 x 4.67 x 4 in.
Power supply unit	12 V, 2 A external adapter
Power supply	100-230 V, 0.3 A, 50-60 Hz
Additional information	2.4" TFT touchscreen



DENTAMILE CURE MC

Equipped with dual-frequency Omni 360° technology, an advanced UV LED arrangement, the DentaMile Cure MC ensures uniform curing through two wavelengths for a variety of printing resins.

TECHNICAL DATA

Capacity	5.70 x 5.70 x 4.13 in.
Heating	Eco-passive up to 60 °C
Interface	2.4" TFT touchscreen
Connection	2.4 GHz WiFi
Power supply	100-230 V, 1.2-2 A, 50-60 Hz
Additional information	<ul style="list-style-type: none"> • Color touch display • Preset exposure programs • LED technology (365 & 405 nm)