

01:00

One-Minute **Product** Meeting

DentaMile Desk RS System

Discover the complete and integrated 3D digital printing solution that enables a more simplified, validated digital workflow.

DentaMile Desk RS-5

The DentaMile Desk RS-5 is a compact, high-precision, high-speed 3D printer designed for a wide range of dental applications using advanced DLP technology. Now enhanced with **RS Turbo** mode, it reduces print times by an additional 20%, while the innovative **Crystal Polish Reservoir** minimizes manual polishing by 90%. The integrated **heating chamber** optimizes print conditions for exceptional accuracy and consistency. With rapid print speeds and **RFID technology** for automatic material identification, the DentaMile Desk RS-5 sets a new standard for efficiency and performance.

DentaMile Wash RS

The DentaMile Wash RS is a fully integrated unit that washes 3D printed devices using a DLP process. The pre-set programs systematically wash the remaining non-polymerized resin from the dental application.

DentaMile Cure RS

The DentaMile Cure RS is a fully integrated curing unit that cures 3D resins with 360° illumination. The pre-set programs lock in the printed mechanical properties ensuring proper curing according to individual resin specifications.

FEATURES & BENEFITS

- Unprecedented quality
- Increased production
- Reproducible precision
- Coordinated processes
- Intuitive operability
- Predictable outcomes
- Guaranteed regulatory compliance

COMPANION PRODUCTS

- LuxaPrint resins
- DentaMile connect

INDICATIONS

- Splints
- Crowns
- Nightguards
- Dentures
- Bleaching trays
- Indirect bonding trays
- Surgical guides
- Impression trays

ORDERING INFORMATION

Item #	Description
170451	DentaMile Desk RS-5 (accessories included)
170452	DentaMile Wash RS
170453	DentaMile Cure RS
170319	3Dmax/DentaMile Desk RS-5 Resin Reservoir
170320	3Dmax/DentaMile Desk RS-5 Resin Reservoir Cover
170321	3Dmax/DentaMile Desk RS-5 Reservoir Storage Cabinet



For more product information visit:
dmg-america.com/products

t: 800.662.6383
e: info@dmg-america.com
w: www.dmg-america.com

A smile ahead together

