

00:01:00

ONE MINUTE PRODUCT MEETING

3Demax, 3Decure, 3Dewash

Discover the state-of-the-art, complete and integrated 3D digital printing solution that enables a more simplified, validated digital workflow.

3Demax

The 3Demax is a compact, high-precision, high-speed 3D printer that can be used for a wide range of applications using digital light processing (DLP) technology. The patented “Force Feedback” technology accelerates printing by up to 50% and RFID technology makes for easy calibration and automatic material identification.

3Dewash

The DMG 3Dewash is a fully integrated unit that washes 3D printed devices using a DLP process. The pre-set programs systematically wash the remaining non-polymerized resin from the dental application.

3Decure

The DMG 3Decure is a fully integrated curing unit that cures 3D resins with 360° illumination. The pre-set programs lock in the printed mechanical properties ensuring proper curing according to individual resin specifications.

FEATURES & BENEFITS

- › Unprecedented quality
- › Increased production
- › Reproducible precision
- › Coordinated processes
- › Intuitive operability
- › Predictable outcomes
- › Guaranteed regulatory compliance

COMPANION PRODUCTS

- › LuxaPrint resins

ITEM # & DESCRIPTION

- 171355 | 3Demax (accessories included)
- 171317 | 3Dewash
- 171316 | 3Decure
- 170319 | 3Demax Resin Reservoir
- 170320 | 3Demax Resin Reservoir Cover
- 170321 | 3Demax Resin Reservoir Storage Cabinet



LOOKING FOR MORE INFO?

Visit www.dmg-america.com/3dprinting