



# Honigum

## The complete impression material system

Trusted by some of the world's top clinicians, Honigum VPS impression material is the ultimate in flowability, hydrophilicity and accuracy. Perfect for crown and bridge, inlay/onlay, and full mouth restorations, its true thixotropic formula does not slump or run, but remains stable where you place it.



# Perfectly Precise Impressions

Honigum is a complete innovative, premium VPS impression material system that offers unmatched detail and exceptional strength to resist tearing upon removal, creating perfect fitting restorations of all types with no need for retakes.

With superior flowability and outstanding stability under pressure, it reaches even the tiniest of crevices to deliver maximum precision impressions. DMG's rapid "snap-set" technology drastically reduces impression errors and distortions due to micro-movement during setting. Plus, it combines a conveniently long working time with a shorter intraoral setting time, resulting in greater control and optimal patient comfort.

**When speed matters**, Honigum's Quad Fast line of materials offer a faster-setting option without compromising superior precision, stability and flawless impression detail.

**For the edentulous application**, Honigum Rigid X-tra Fast is a fast-setting, rigid, heavy bodied VPS impression material for customized impression trays, including edentulous impression techniques. Easily removed, Mono Quad Fast works well with Rigid X-tra Fast as a wash material, providing superior accuracy and greater muscular detail without compressing soft tissue.


## Ordering Information

### Honigum

Item#	Product	Contents
999539	Honigum Light	4-25ml Cartridges, 16 Automix and Intra-Oral Tips
999831	Honigum Light	2-50ml Cartridges, 10 Automix and Intra-Oral Tips
999639	Honigum Light Quad Fast	4-25ml Cartridges, 16 Automix and Intra-Oral Tips
999832	Honigum Light Quad Fast	2-50ml Cartridges, 10 Automix and Intra-Oral Tips
999830	Honigum Mono	4-50ml Cartridges, 12 Automix and Intra-Oral Tips
999839	Honigum Mono Quad Fast	4-50ml Cartridges, 8 Automix and Intra-Oral Tips
999766	Honigum Heavy	4-50ml Cartridges, 12 Automix Tips
999838	Honigum Heavy Quad Fast	4-50ml Cartridges, 12 Automix Tips
999900	Honigum Rigid X-tra Fast	4-50ml Cartridges, 8 Mixing Tips
999568	Honigum Mono	1-380ml Cartridge, 10 MixStar Tips
999537	Honigum Heavy	1-380ml Cartridge, 10 MixStar Tips
999638	Honigum Heavy Quad Fast	1-380ml Cartridge, 10 MixStar Tips

### Tips and Accessories

Item#	Product	Contents
999482	MixStar eMotion Machine	1-Mixing Machine Unit
999567	MixStar White Mixing Tips	Box of 50 Tips
999400	Honigum Automix Tips	Bag of 50 Tips (Heavy, Mono)
110424	Automix/Intra-Oral Mixing Tips	Bag of 35 Each (Light)
999490	Intra-Oral Tips	Bag of 50 Tips (Mono)
999394	Tray Adhesive	1-10ml Bottle
110253	Type 25 1:1 Applicator Gun	1-Gun
999507	Type 50 1:1 Applicator Gun	1-Gun



DMG America, LLC  
 65 Challenger Road, Suite 340, Ridgefield Park, NJ 07660  
 tel: 800.662.6383 email: info@dmg-america.com  
 web: www.dmg-america.com

