



Safety-Wipes™

Sanitary. Smart. Secure.

A smile ahead together



Safety-Wipes™

Elevate Your Clinical Workflow

Whether you're in hygiene, restorative, or endo, Safety-Wipes™ help you streamline every procedure with added safety and precision.



Developed to reduce hand-stress and facilitate procedures when unassisted, these adhesive-backed foam pads stick where you need them—gloves, bibs, or trays—so you can focus on patient care, not fumbling with gauze or stopping mid-procedure due to punctures.

Safety-Wipes Features & Benefits

Hands-Free Convenience

Sticks directly to your glove or tray—always within reach

Superior Protection

Thick, dense foam reduces puncture risk and protects against sharps injuries

Time-Saving Efficiency

No more reaching for gauze mid-procedure—wipe, clean, or blot with one hand

Multi-Functional Design

Use to hold files, absorb fluids, or apply antiseptic in canals or grooves

Ergonomic Relief

Reduces finger strain from squeezing or holding traditional gauze pads

Ordering Information

Item #	Product	Contents
163020	Safety Wipes	100 wipes per bag (2 bags of 50)

DMG America, LLC
65 Challenger Road, Suite 340, Ridgefield Park, NJ 07660
tel: 800.662.6383 email: info@dmg-america.com
web: www.dmg-america.com

A smile ahead together

